# Norav medical Inc.

---

CE 0473

## PC-ECG 1200

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


## ZUBEHÖR INSTALLATION

D/A Board: Model: CIO-DAC02 for ISA bus

JUMPERS SETTING: D/A0= -5, D/A1= -5, WAIT STATE= OFF

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Board P/N and Board target	board address and address switches	COMMON pin # on D25 connector	OUT 1 pin # on D25 connector	OUT 2 PIN # on D25 connector
D1-a Treadmill Analog control	300H. switches 8,9 DOWN the rest UP	1	Speed pin : 24	Grade pin: 18
D1-a Ergometer Analog control	300H. switches 8,9 DOWN the rest UP	1	Power pin : 24	
D1-s Metabolic Interface	304H. switches 8,9,2 DOWN the rest UP	1	Heart Rate pin: 24	Work Load pin : 18
D1-t ECG TRIGGER	308H. switches 8,9,3 DOWN the rest UP	1	TTL OUT pin: 24	

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Board P/N and Board target	board address and address switches	COMMON pin # on D25 connector	OUT 1 pin # on D25 connector	OUT 2 PIN # on D25 connector
D1-a Treadmill Analog control	300H. switches 8,9 DOWN the rest UP	1	Speed pin : 23	Grade pin: 17
D1-a Ergometer Analog control	300H. switches 8,9 DOWN the rest UP	1	Power pin : 23	
D1-s Metabolic Interface	304H. switches 8,9,2 DOWN the rest UP	1	Heart Rate pin: 23	Work Load pin : 17

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Connector pin 21 an D37 connector – common

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**3X4 FORMAT:** Das klassische Format. 12 Kanal EKG von 2.5sek EKG + 10sec Streifen that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

 **6X2 FORMAT:** 12 Kanal EKG von 5sec EKG + 10sek Streifen. Note that according to SE the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

 **12X1 FORMAT:** 12 lead ECG of 10sec ECG. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead

 **3X1 FORMAT:** 3 lead ECG of 10sec ECG. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

**Mittlung:** Produces a typical QRS for every lead from the raw ECG data.

The averages can be printed either full size or minimized in a credit card size. To perform miniature print enters FILE and that PRINT MINIATURE.

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**50/60 Hz:** AN/AUS for line interference filter. The user should OPTIONS for 50 or 60 Hz prior to operation.

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- ? Record a rest study and apply MATLAB function under FILE menue.
- ? A file with extension LIB will be created in 'default data folder'.  
The file has the following structure:

10,000 bytes X 12 Leads (I,II,III,AVR,AVL,AVF,V1,V2,V3,V4,V5,V6);  
1sec = 1000 bytes = 500 samples.

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For example if lead AVR (its number is 4 ) from the second 3 is requested:

`fseek(handle,3*10000+2*500,-1)`  
`t = fread(handle, [sec * 500 ,1], 'int16')`  
where sec(1,2,3,4,...,10)

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### Explanation on the Metabolic Stress Estimation:

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to be used when PC-ECG 1200 shares the same PC with another application in real time

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Hersteller	Modell		RS232 connector	RS232 type**
Trackmaster	TM425		D-9 male	Straight
Trackmaster	TMX425		D-9 male	Straight
Parker	PM		D-9 female	Straight
RAM	770		D-9 female	Straight
Woodway	PPS55-MED		D-9 male	Crossed
Powerjog	GM, J			
Biomedical Sys				
HP COSMOS	MERCURY-MED			
QUINTON	TM-55			
QUINTON	ST-55			*

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seitig und pin 5 PC seitig verbindet pin 5 Laufband seitig. . Crossed type means that pin 2 on PC side connects to pin 3 on treadmill side, pin 3 on PC side connects to pin 2 on treadmill side and pin 5 on PC side connects to pin 5 on treadmill side.

? In order to determine if a certain treadmill requires a straight or crossed cable one can see on the treadmill connector what pins are assigned for TD and RD. If pin 2 is RD and pin 3 is TD it requires a crossed cable if vice versa it requires a straight one.

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 **Bemerkungen:** Enter free text during or after the ECG recording. Will be printed and saved out together with ECG traces. 

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<b>Siehe unter:</b>	<b>Funktion</b>	<b>Beschreibung</b>
<b>EKG Aufzeichnung</b>		
	Speicher Optionen	If "Auto Save" is ON the file will be stored by Last name or by ID. If "Auto Save" is OFF the program will ask if to store and the user will be asked to enter filename prior to storing.
	Simulator EKG	If unmarked (default) ECG recording is done from PC-ECG unit. If marked the ECG recording is done from demo file included in the software package. In this case the recording unit is not needed.
	Stop to confirm QRS	If ON the user can choose the "normal QRS" else the program will make automatic decision.
	Template Correlation	Allows to define what percentage the QRS matching during signal averaging will be. The higher the number the better the matching is
	Target Number	Number of typical heartbeats that will be counted during the averaging stage.
	Data Directory	Enables the user to define his own directory for saved ECG recordings
<b>Zeige</b>		
	ECG Colors	Choose preferable colors for monitoring window.
	Average Colors	Choose preferable colors for averages and filtered ECG windows.
	HR Trend Colors	Choose preferable colors for Heart Rate window.
	Restore Defaults	When activated the Factory defaults are restored.
<b>Installation</b>		Saves users data of Hospital and Physician. This data will be printed on any print out and send as mail.
<b>Ausrüstung</b>	COM PCECG	Determines the serial input for the PC-ECG unit
	Graph paper	When ON it prints 1mm and 5 mm squares on printouts. Regular grid is guaranteed to fit any printer. Improved grid shows a fine grid but may not work on some printers.
	Display Size	Choose between 14/15-inch screen (default) and 17 inch screen. This setting is required in order to display the ECG and grid in correct scale.

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**START/STOP MITTLUNG.** Use this button to start the averaging period. The averaging period default is 200 beats. It can be changed in OPTIONS, ECG RECORDING, and TARGET NUMBER.

**GITTERNETZ:** optional display 5mm raster. Print outs are always with 1mm raster.

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Upon completion of the signal-averaging phase the result screen appears:

On the upper side from left to right are displayed: signal averages, filtered averages, Vector Magnitude (on which the calculations are made). On the lower hand side are displayed numerical results and Heart trends.

The numerical results are:

- ? Gemittelte Schläge: Anzahl der gemittelten und bewerteten Normalschläge, die während des Tests aufgezeichnet wurden.
- ? HF Mittlung in Schläge/min
- ? Hochfrequenz QRS Duration in milliseconds from ONSET to OFFSET. This is the width of the filtered QRS containing only high frequencies. The higher this numbers the higher the patient risk.
- ? RMS LASTS 40 milliseconds in Mikrovolt. This parameter informs what is the total activity at the last portion of the QRS. The lower this numbers the higher the patient risk.
- ? DURATION UNTER 40 Mikrovolt in Millisekunden. Dieser Parameter zeigt, wie stiff the offset of the QRS is. Je höher diese Anzahl, desto höher ist das Risiko des Patienten.
- ? H.F. Rauschen: in microvolts. This information is for quality assurance of the study. The lower these number the higher the result accuracy.  
The results are calculated automatically. The user can overrule the automatic positioning of the ONSET/OFFSET markers by means of direction keys.

## Monitoring

This option package enables long term recording and storage into the disk. The user can decide what leads and in which sample rate he likes to monitor on the screen and save into the disk. During the study real time printing on thermal printer is possible.

### Icons:

**NEW RECORDING: Start** new monitoring. The patient data can be entered prior to ECG recording, but it is not mandatory. The recording time goes according to the SETUP for ECG RECORDING. The user can stop recording by hitting the GO/STOP icon.

**OPEN RECORDING:** Shows recordings that were saved on disk.

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**PATIENTDATEN:** Patientendaten Information.

**BEMERKUNGEN:** Enter free text during or after the ECG recording. Will be printed and saved out together with ECG traces.

**SETUP:** Enable the customer to tailor his own preferences for operation.

<u>See under:</u>	<u>Function</u>	<u>Description</u>
<b>Record Leads Setting</b>	Leads	Choose requested leads.
	Sample rate	Choose requested samples per second per channel
<b>12 Leads View</b>	Default triplet	This function is operational only if 12 leads monitoring was chosen. 3 leads to appear if 3X1 format is used
	Strip Lead	10 sec lead to appear at 4X3 and 6X2 formats
<b>ECG recording</b>	Filter 50/60Hz	When marked the default status of 50/60Hz filter will be ON (according to the marked frequency 50 or 60). Default is unmarked.
	EMG Filter	When marked the default status of the EMG filter is ON. Default is unmarked.
	Save options	If "Auto Save" is ON the file will be stored by Last name or by ID. If "Auto Save" is OFF the program will ask if to store and the user will be asked to enter filename prior to storing.

	Auto stop nach 10 sek	If unmarked (default) recording run till stopped by user if marked stops recording automatically after 10 sec.
	Simulator EKG	If unmarked (default) ECG recording is done from PC-ECG unit. If marked the ECG recording is done from demo file included in the software package. In this case the recording unit is not needed.
	Ekg Aufnahmezeit (h:m)	Determines study duration in Minuten.
	DatenDirectory	Enables the user to define his own directory for saved ECG recordings. In case hard drive D: exist it is advisable to us it in order not to consume drive C: space.
<b>ZEIGE</b>	Draw over lead borders	When marked (default) it does not limit the extreme high amplitude ECG pulses to exceed the borders. When unmarked it will chop the pulses at the borders.
	Horizontal Scale	Sets the default value for the Horizontal scale window on the screen.
	Vertical Scale	Sets the default value for the Vertical scale window on the screen.
	Slider step size	Off line function. Sets the default value for the slider when moved by mouse or arrow keys.
	colors	Enables the user to choose his own colors.
	Restore Defaults	When activated the Factory defaults are restored.
<b>Installation</b>		Saves users data of Hospital and Physician. This data will be printed on any print out and send as mail.
<b>Ausrüstung</b>	COM PCECG	Determines the serial input for the PC-ECG unit
	Display Size	Choose between 14/15-inch screen (default) and 17 inch screen. This setting is required in order to display the ECG and grid in correct scale.
	Graph paper	When ON it prints 1mm and 5 mm squares on printouts. Regular grid is guaranteed to fit any printer. Improved grid shows a fine grid but may not work on some printers.
	Large Remarks font	Allows printed text to appear larger
	Farbausdruck	If color printer is used it allows forcing the printer to print B/W only.
	Thermo Plotter	Sets LPT port for optional thermal paper

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**3X4 FORMAT:** The classical format. 12 lead ECG of 2.5sec ECG + 10sec trace. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

**6X2 FORMAT:** 12 lead ECG of 5sec ECG + 10sec trace. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

**12X1 FORMAT:** 12 lead ECG of 10sec ECG. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead

**3X1 FORMAT:** 3 lead ECG of 10sec ECG. Note that according to SETUP the ECG data can be limited in amplitude or unlimited which can cause to one lead data to overlap a neighboring lead.

**NEXT LEADS:** Scroll for all leads in the 3X1 display

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**50/60 Hz:** ON/OFF for line interference filter. The user should OPTIONS for 50 or 60 Hz prior to operation.

**EMG:** ON/OFF für Muskelrauschfilter.

HORIZONTAL SCALE: Shows scale in mm/sec

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Cur. Time: Shows the time where the slider is positioned

Testzeit: Entire study time

### Operation by Functional Keys:

**F1- New Recording**

**F2- Start/Stop**

**F6 – Print**

**F11- Open saved Study**

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**Under FILE**

**Export to rest - A piece of 10 sec containing original leads I,II,V1-V6 and calculated leads III, aVR, aVL, aVF is transferred into rest format (up to 12 leads 10 sec). Calculated leads are performed only if I and II are acquired.**

**Export to MATLAB - A piece of 10 sec containing acquired leads is transferred into MATLAB format .**

**Import from ISHNE – Long term high resolution ECG recorded on Holter can be transferred into a monitoring study.**

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Time and Frequency Domain analysis designed for Short Studies in which one or several time segments are measured ; as in a Tilt Study.

Measured/Reported parameters according to NASPE/ESC Guidelines.

Upon hitting the HRV icon the working screen opens. It is recommended to review the default set up. Hit the hammer icon to get the following folders:

Sample Rate/ Leads Selection – choose leads and resolution.

ECG Recording – Set filters to be active, define if ECG should be saved during the study, set for target number of beats. If the data base should be used check “Use ECG Database” otherwise define where to save studies by using “Browse” if a different location than the default is needed.

View – If default colors need to be changed.

Installation – Enter user details and change to USA for USA height and weight.

Environment – Enables to set according to used PC.

Start the study by hitting “new” and entering patient data. The ECG leads are monitored on the screen for quality check. If satisfied, hit the R-R icon. The ECG leads are now shown on the upper part of the screen and a Tachogram trends built up in the middle strip. On the bottom a slider shows the study status and time. During the study a new time segment (interval) can be defined by hitting the flag icon (interval). Each interval must be named during the study in order to remain as a valid interval.

When all predefined beats are done or if terminated by hitting the R-R icon the following screen appears:

Upper left hand side: Histogram related to the active part (yellow) of the tachogram

Upper right hand side: Power Spectrum Distribution

Middle Left hand: The tachogram trend showing all intervals. Each beginning of interval is marked with a red line followed by the interval’s name. To activate an interval click on it with the mouse. In order to activate several neighboring segments use the slider while “Shift” key is pressed.

The slider on the bottom is used to define new intervals, change interval duration and activate several intervals.

In order to add or substrate an interval use the + and – on the middle right side in conjunction with the slider. The pencil icon is used to edit interval’s names.

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max RR - längste R-R Periode

min RR - kürzeste R-R Periode

Average RR – Average of interval in active interval

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RMSSD – Die Quadratwurzel of the mean of the sum of the squares of differences between adjacent NN intervals.

HRV triangular Index - Total number of all R-R intervals divided by the height of the histogram of all R-R intervals measured on a discrete scale.

### Frequency Domain:

ULF – Power of the ultra low frequency range

VLF – Power of the very low frequency range

LF – Power of the low frequency range

HF – Power of the high frequency range

### **HFV Intervall Messung:**

A QRS detector measures interval between any 2 valid beats. It calculates a sliding N-N average and compares each interval to it. When a significant change occurs, the current beat is either a premature (as in PVC) or a prolonged which may indicate either a compensatory pause following a PVC or a missing beat. A premature interval and a following prolonged interval (compensatory pause) timed in the range of twice the current N-N interval are averaged. This methodology maintains a consistent time axis in the presence of PVCs. If a prolonged interval follows a normal interval, but at twice the current N-N interval, it indicates the presence of a missing beat. Accordingly, the missing beat is computed as present. The last event is very rare, because the recording is done in rest condition and the software has no problem to detect the QRS in such condition

## **EKG DATENBANK**

This option enables to keep and manage ECG studies in a catalog which is oriented according to patient name or ID. In case of network all users can share the same data base (which is recommended to locate in the server). All data acquisition for all applications can be initiated either in the application itself (as previous versions) or from the data base main screen.

Upon entering to the data base for the first time the user will be prompt to confirm the location of the data base. Choose NEW if no data base exist already and OPEN if one exist.

The screen shows 2 areas: On the left hand side a list of patients that consist of 3 columns of ID , Last Name, and First Name. The list can be ordered according to alphabetical order of each column. On the right hand side is the studies information area. When a certain patient is marked in the list by means of selecting by Find or by mouse the studies contained for this patient will be displayed.

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Import – Add studies (windows and DOS types) that were recorded with out the usage of the data base. In order to select all files of a directory which contains patient data type CTRL A

Setup – Enables to define the location of default Workspace, Patient identification, and a special file called External File. This file in a format of windows.INI files allows to prepare a list of patients that can be read by the PC-ECG 1200 applications. Note : this feature does not require the D1 software key.

About – Keeps version number, Communication and web information and disk and RAM memory size. The HASP ID number is used for identification of the software key for the purpose of adding additional software options into the key.

Gruppen – Is used to define different groups that may be used to relate a patient to. Wie Privat, HMO, Militär etc.

Suche – Patient can be find by entering a string into any or all fields of ID , Last Name and First Name.

Öffne – Before using a patient must be marked. The patient information is shown.

Neu – Insert new patient. The User is prompt to enter part or all details of the patient. It is recommended to enter at least ID, Last and First Name. If patient details matches an existing one the addition of this patient to the list will be prohibited. The patient that matches the details will be marked to allow the user to add a study to that patient.

Lösche – Ability to delete an entry. If the entry is not empty of studies the user will be asked to confirm.

Test – Go to application and start recording.

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Property – Shows study properties

Check – Is used to select and show on the screen only studies that have been marked in the left hand box.

Kopie – Ability to copy a study.

Move – Moves a study and deletes it from data base. The default option is to leave the study's property in the data base. The user can decide to remove the entire study. If the default option was used then upon completion of the operation the study will remain in the list with a X sign.

Lösche - Deletes a study from data base. The default option is to leave the study's property in the data base. The user can decide to delete the entire study. If the default option was used then upon completion of the operation the study will remain in the list with a X sign.

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<b>Merkmal</b>	<b>Modell</b>		
	<b>1200B</b>	<b>1200S</b>	<b>1200M</b>
Grösse [cm]	15x12x2	20x14x3.5	15x8x3
Gewicht [gr]	250 (inclusive Gurt/Tasche)	500	200
EKG samples/sek	250, 500 ,1000,2000		
A/D bits	12 ( 2.44 µV/LSB)		
Defibrillationsschutz	eingebaut		Mit Kabel C1-d
Simultan 12 Kanal	alle		
CMMR	> 100 dB		
Eingangs Impedance	> 100Mohm		
Signal dynamic range	10mV		
DC max. Eingang	+ oder - 330mV		
Frequenzbreite (-3db)	0.05 - 300 Hz		
Low Pass Filter (Software)	35 Hz		
Base Line (software)			
Line noise Filter (software)	50/60Hz		
Sicherheitsstandards	IEC 601-1 , IEC 601-2-25 , IEC 601-2-27		

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Zeige Programmversion und copyright